

Turn on by pressing control knob. Click and turn up level a little before doing 3 button trick



- 3 BUTTON TRICK**
1. Press and hold control knob
 2. Press pause and test buttons
 3. Let all three go.

Diagram 1

Start
Put the ODFS Pace into SETUP MODE using three button trick

Find the head of fibula
If skin is dry, moisten with water
If skin is greasy clean with an alcohol wipe

Place electrodes in standard position
Patient should be seated with leg supported with and extended knee

Select NEW SETUP?
Select DROPPED FOOT
Select HEEL RISE
Set the OUTPUT CURRENT



Increase the current until dorsiflexion with eversion is seen. The pulse width is set automatically to 50%

Correct movement Dorsiflexion with eversion?

Sensation OK?

Correct movement Dorsiflexion with eversion?

If you nearly have sufficient dorsiflexion and eversion try walking anyway. Often the level can be increase when standing as sensation is often reduced



Anatomy

Common Peroneal Nerve
Tibial nerve
Superficial branch
Eversion and Plantarflexion
peroneus longus
peroneus brevis
Deep branch
Dorsiflexion and Inversion
tibialis anterior
extensor digitorum longus
extensor hallucis longus
peroneus tertius


Press and release the test button. Repeat as required



The current can only be turned up while the LED flashes. Can turn down at any time.

Use exercise stimulation. building up level day by day until acceptable 2x20min a day for 1 month

Popliteal Fossa position for the withdrawal reflex (knee flexion) and a stronger effect



Electrode on the Biceps femoris tendon and lateral boarder of the popliteal fossa
Avoid the tibial nerve

Too much inversion?

Move the active electrode (black lead) backwards and upwards about 1cm

Still too much inversion?

Move active electrode (black plug) to popliteal fossa

Too much eversion?

Move the active electrode (black lead) forwards and downwards about 1cm

Still too much eversion?

Swap the active (black plug) and indifferent (red plug)

Eversion but no dorsiflexion?

Increase stimulation level

Weak response?

Increase current

Increase OUTPUT FREQ

8th menu option OUTPUT FREQ.

Increase in 50ms steps



Spastic calf reduces ROM or causes clonus spasm?

Increase rising edge ramp NB ramp must be quick enough to pick up the foot in walking


2nd menu option RISING RAMP

Increase in 50ms steps



Reduced knee flexion, calf co-contraction

Reduce Stimulation frequency



Dorsiflexion without eversion
Motor point of anterior tibialis
Avoid tibia bone
Motor point stimulation requires high levels of stimulation and may be more uncomfortable

Put active (black plug) on to the motor point of Tibialis Anterior

Still too much eversion?

To much inversion?

Change waveform to Symmetrical Biphasic

Correct movement Dorsiflexion with eversion?

Repeat analysis and try again

If the correct movement is still not achieved the ODFS Pace may not be a suitable treatment for this Patient

Prepare for walking
Go to diagram 2 (Stay in set up mode)

7th menu option OUTPUT WAVEFORM

