

The Odstock Medical Ltd Pace Sleeve is designed to cover and hold your stimulation electrodes in place while also providing a pocket to hold the ODFS® Pace (XL) stimulator on your leg.

The Pace Sleeve, when set up on your leg, enables the use of shorter electrode leads (25 cm) and footswitch leads (35 cm) for convenience and efficiency.

The Pace Sleeve is supplied in packs of 2.

Code	Size	Approximate range of leg circumference
01-005-0049	L130mm x W80mm	200mm to 260mm
01-005-0050	L150mm x W100mm	250mm to 350mm
01-005-0051	L175mm x W120mm	330mm to 430mm
01-005-0052	L175mm x W140mm	430mm to 500mm
01-005-0053	L175mm x W160mm	460mm to 560mm



\*while stocks last



25cm electrode leads and 35cm footswitch leads are available for wearing the ODFS Pace on the leg

To order, please contact [sales@odstockmedical.com](mailto:sales@odstockmedical.com) or call 01722 439540

11-005-0022 v3

Email: [enquiries@odstockmedical.com](mailto:enquiries@odstockmedical.com)

Call: 01722 439 540

Website: [www.odstockmedical.com](http://www.odstockmedical.com)

Videos: [Odstock Medical Ltd](#)

Follow us: @OdstockMedicallimited @OdstockMedical



Scan here