

Odstock Medical Limited (OML) produces innovative FES device upgrades and accessories to enhance FES treatment. Some of these options may not be covered under NHS funded treatment so may require personal payment to obtain them.



ODFS® Pace XL

The ODFS® Pace XL uses a wireless connection between the OML LINQ™ footswitch and the stimulator. This can make the device easier and more convenient to use. Existing ODFS® Pace devices can be upgraded to an ODFS® Pace XL (only available to patients of Salisbury, Outreach and UCL who have devices issued after September 2021). Please contact OML to find out more.



OML LINQ™

The OML LINQ™ is a wireless transmitter that relays the footswitch signal to the stimulator. This removes the need for the connection cable. The OML LINQ™ requires a coin cell battery, which needs replacing approximately every three weeks.



If the current FES treatment is being funded by a local funding body, upgrading would require agreement from them or alternatively the user can pay for it. If FES treatment is being self-funded, the user would need to fund the upgrade. The upgrade fee is a one-off charge and also requires the purchase of the accessory.



OML Remote

The OML Remote is a hand-held remote control for the ODFS® Pace XL, allowing the stimulator to be easily and conveniently controlled. The OML Remote can be used to start and stop stimulation, change the intensity, test the stimulation, switch between walking and exercise modes and check the battery status.

The Remote also allows two ODFS® Pace XL stimulators to be controlled with a single OML LINQ™ footswitch, enabling applications such as correction of bilateral dropped foot and correction of dropped foot with hip, knee or calf control.

What our FES users say...

“ It's impossible for me to overstate the benefits FES has for my life. ”

11-005-0030 v3.0

Email: enquiries@odstockmedical.com

Call: 01722 439 540

Website: www.odstockmedical.com

Videos:  [Odstock Medical Ltd](#)

Follow us:  @OdstockMedicalLimited  @OdstockMedical



Scan here