



OML

Leading Rehabilitation
Through Technology

FES User Day 15th June 2018



OML LINQ™ Refresher

Stacey Finn





“The LINQ has enabled me to wear any shoes”

“Much easier to use than wired footswitch for a hemiplegic, much more easier and comfortable than the wireless insole”

“For me the LINQ is essential for wearing a dress, skirt and shorter trousers”

- Easy to use
- Commonly clipped to side of shoe or tucked in sock
- Commonly moved to different shoes
- Used with various shoe types
- Users would recommend to other FES users

- Some challenges joining. Button intentionally hard to press. Don't forget quick-start and videos on web-site.
- Do not need to re-join to Pace XL when changing battery.
- Hook one side of the LINQ into the clip first.
- 2 Footswitch sockets available with the application of running and stair climbing.
- Rubber coating can mark with abrasion.